

AstraZeneca selected monitors

- Solution is contactless and driverless (platform independent)
- Supports wireless and USB loggers

Sentinel



- AZ approved temperature monitor to be used for the release of AZ products
- Re-usable monitor
- Wireless connection (Bluetooth) to Sentry gateway

Sentry600 Gateway



- Fixed Gateway to be permanently connected to a power outlet.
- Transfers temperature data to CCIS over mobile network

Sentry600 In-Transit



- Re-usable In-Transit near real-time GPS monitor
- Connects with Sentinel during transit to transmit temperature and location data to CCIS over mobile network

webLogger II



- Single use temperature monitor
- USB connection to computer

webLogger II dry ice



- Single use temperature monitor for deep frozen
- USB connection to computer